**ANTI-NEUTROPHIL CYTOPLASMIC ANTIBODY (P-ANCA)
TEST SYSTEM**

For *In Vitro* Diagnostic Use

AD PAN60 60 Tests

INTRODUCTION:

Standard IFA methods allow for the observation of several different patterns. Two patterns that have been well defined are C-ANCA and P-ANCA. The P-ANCA pattern shows an uneven granular, staining of the cytoplasm. During the 2nd International ANCA Workshop it was agreed that these two different patterns should be used to subclassify the antibodies. ANCA antibody detection by IFA methods has been a useful aid in the assessment of patient diagnosis, and to a certain extent their prognosis and response to therapy. Although there has been much work done with the IFA ANCA, the availability of an EIA MPO Test has further enhanced the overall picture. The use of EIA and IFA together allows for the best possible patient assessment. The MPO & ANCA (PR3) EIA allows for a very quick qualitative as well as quantitative report.

PRINCIPLES:

The primary reaction in this assay involves human antibody (patient sera) and a specific antigen (human granulocytes). If ANCA antibody is present in the patient sera it will bind to form an antigen/ antibody complex. This complex is then labeled with an FITC labeled antihuman conjugate that allows one to visualize the reaction through the microscope.

MATERIALS PROVIDED:

1. FITC IgG H&L conjugate No. AD CGEA25 (2.5 ml) w/Evans Blue is to be stored at 2-8° C upon receipt. The conjugate is stable at this temperature until expiration date on label.
2. The Formalin Fixed antigen slides of ANCA (Human Granulocyte) substrate **must be stored at 2-8° C** upon receipt. Check label for specific expiration date.
3. P-ANCA positive control No. AD PCPA (1.0 ml), demonstrating a P-ANCA pattern, should be stored at 2-8° C upon receipt. Check label for specific expiration date.
4. ANCA negative control No. AD NC (1.0 ml) should be stored at 2-8° C upon receipt. Check label for specific expiration date.
5. Buffer Pack No. AD PBS1 - Phosphate Buffered Saline is stable at room temperature storage. Check label for specific expiration date. Check label for specific expiration date. Rehydrate buffer with 1 liter of sterile DI H₂O. The reconstituted, buffer does not contain preservatives and should be stored at 2-8° C. Care should be taken to avoid contamination.
6. Mounting Medium No. AD AMM3 should be stored at 2-8° C. Check label for specific expiration date.

GENERAL STORAGE OF THE KIT:

The slides, controls, and conjugate should be stored at 2-8° C.

ADDITIONAL MATERIALS REQUIRED BUT NOT PROVIDED:

- Test tube rack or microtiter system
- Disposable pipettes
- Staining Dish and Slide Forceps
- Fluorescence Microscope
- Moisture Chamber
- Volumetric Flask (500 ml)
- Distilled H₂O

ADDITIONAL COMPONENTS AVAILABLE:

C-ANCA Positive Control (Cat. AD PCCA)

Note: All kit components are available separately. Please see catalog for more details.

SPECIMEN COLLECTION:

Serological specimens should be collected under aseptic conditions. Hemolysis is avoided through prompt separation of the serum from the clot. Serum should be stored at 2-8° C if it is to be analyzed within a few days. Serum may be held for 3 to 6 months by storage at -20° C or lower. Lipemic and strongly hemolytic serum should be avoided. When specimens are shipped at ambient temperatures, addition of a preservative such as 0.01% thimerosal or < 0.1% sodium azide is strongly recommended.

TEST INSTRUCTIONS:

Screening: Dilute test serums 1:20 (1 part patient sample to 19 part diluent) in PBS.

Titration: Set up doubling dilutions of serum starting at 1:20, (i.e. 1:20, 1:40, 1:80, 1:160, 1:320, etc.). The slides are ready to use after they reach room temperature.

1. Allow slide to reach room temperature before opening envelope. Tear envelope at notch. Carefully remove the slide and avoid touching the antigen areas. The slide is now ready to use.
2. Place a drop of serum (20-30 µl) over the antigen wells.
3. Place slide with patient's serum and controls in a moist chamber for 20 minutes at room temperature (approximately 19-24° C).
4. Remove slide from moisture chamber and allow the serum to run off onto a piece of paper towel. Using a wash bottle, gently rinse remaining sera from slide being careful not to aim the stream directly on the well.
5. Wash in PBS (Cat. 1601) for two separate five minute changes.
6. Place a blotter on the lab table with absorbent side up. Remove slide from PBS and invert so that substrate side faces absorbent side of blotter. Line up wells to blotter holes. Place slide on top of blotter. Do not allow substrate to dry.
7. Wipe back of slide with dry lint free paper towel. Apply sufficient pressure to slide while wiping to absorb buffer. Do not allow substrate to dry.
8. Deliver 1 drop (20-30 µl) of conjugate per antigen well. Repeat steps 3-7.
9. Place 4-5 drops of mounting medium on slide.
10. Apply a 22 x 70 mm coverslip. Examine the slide under a fluorescent microscope.

Note: To maintain fluorescence, store mounted slide in a moisture chamber placed in a dark refrigerator.

Note: Caution should be taken to not extend incubation or rinse times. The substrate will be effected and poor morphology will result.

RESULTS:

Positive

A positive result is reported when the cytoplasm of the human granulocyte substrate displays a 1 + or greater fluorescence. P-ANCA and C-ANCA will give a similar uneven granular staining of the cytoplasm, with formalin fixation.

C-ANCA and P-ANCA may occur together C-ANCA antibodies are associated with classic Wegener's granulomatosis.

P-ANCA (MPO) antibodies are associated with renal limited diseases.

Negative

A serum is considered negative for ANCA if the cytoplasm fluorescence is less than 1 +. Patients should be screened on ANA HEP2 substrate to avoid confusion with PSEUDO-ANCA. PSEUDO-ANCA will stain the cytoplasm of HEP2 cells whereas true ANCA will be negative on HEP2 unless the patient possesses both ANA and ANCA antibodies.

QUALITY CONTROL:

1. Positive and negative serum controls must be included in each day's testing to confirm reproducibility, sensitivity and specificity of the test procedure.
2. The negative serum control should result in little (+) or no fluorescence. If this control shows bright fluorescence, either the control antigen, conjugate or technique may be at fault.
3. The positive serum control should result in bright 3+ to 4+ fluorescence. If this control shows little or no fluorescence, either the control, antigen, conjugate or technique may be at fault.
4. In addition to positive and negative serum controls, a PBS control should be run to establish that the conjugate is free from nonspecific staining of the antigen substrate. If the antigen shows bright fluorescence in the PBS control repeat using fresh conjugate. If the antigen still fluoresces, either the conjugate or antigen may be at fault.

PRECAUTIONS:

1. All human components have been tested for (HB_sA_g) and HTLVIII/LAV by an FDA approved method and found to be negative (not repeatedly reactive). However, this does not assure the absence of HB_sA_g or HTLVIII/LAV. All human components should be handled with appropriate care.
2. The sodium azide (<0.095%) included in the controls and conjugate is toxic if ingested.
3. Do not use components beyond their expiration date.
4. Follow the procedural instructions exactly as they appear in this insert to insure valid results.
5. For *In vitro* Diagnostic Use
6. Handle slides by the edges since direct pressure on the antigen wells may damage the antigen.
7. Once the procedure has started do not allow the antigen in the wells to dry out. This may result in false negative test results, or unnecessary artifacts.
8. Use separate pipette tips for each sample and reagent to avoid cross contamination.
9. Reagents should be inspected for evidence of bacterial or fungal contamination.
10. Do not reuse substrate slide.

Component	AD PBS1 PBS Powder Packets AD AMM3 Mounting Medium	Precautionary Statement Prevention: P264 Wash thoroughly after handling. P280 Wear protective gloves and clothing.
Pictogram		Response: P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if easy to do. Continue rinsing.
Signal Word	WARNING	
Hazard Statement	H319 Causes serious eye irritation	P337+P313 If irritation persists, get medical advice/attention.

BIBLIOGRAPHY:

1. Wilk A and van der Woude F: Anti-neutrophil cytoplasmic antibodies (ANCA): A historic review. *APMIS* 97 (Suppl 6):7, 1989.
2. Wilk A: Delineation of a standard procedure for indirect immunofluorescence detection of ANCA. *APMIS* 97 (Suppl 6):12-3, 1989.
3. Rasmussen N and Wilk A: Indirect immunofluorescence examination for IgG-ANCA in sera submitted for the 1st international workshop on ANCA, 1988. *APMIS* 97(Suppl 6):16-20, 1989.
4. van der Woude FJ, et al: Autoantibodies against neutrophils and monocytes: Tools for diagnosis and marker of disease activity in Wegener's granulomatosis. *Lancet*, Feb 13:425-9, 1985.
5. Falk R.J. and Jennette J.C.: Anti-neutrophil cytoplasmic autoantibodies with specificity for myeloperoxidase in patients with systemic vasculitis and idiopathic necrotizing and crescentic glomerulonephritis. *N Engl J Med* 38(25):1651-7, 1988.
6. Goldschmeding R., et al: Different immunological specificities and diseases associations of c-ANCA and p-ANCA. *Neth J Med* 36(3):121-5, 1990.

ALPHADIA sa/nv
DIAGNOSTIC PRODUCTS
Avenue Vésale 26
B1300 WAVRE
BELGIUM
TEL : 32 (0) 10 24 26 49
FAX : 32 (0) 10 24 55 99
contact@alphadia.be

EC REP Emergo Europe
Princessegracht 20
2514 AP The Hague
The Netherlands

ANTI-NEUTROPHIL CYTOPLASMIC ANTIBODY (P-ANCA) TEST SYSTEM

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