## ANTI-SKIN ANTIBODY TEST SYSTEM

For *In Vitro* Diagnostic Use

AD SMO48 48 Tests  
AD SMO96 96 Tests

### Introduction:

The *in vitro* detection of skin antibodies by the indirect immunofluorescent technique has been established as an aid in the diagnosis of skin and systemic diseases. The utilization of monkey esophagus, has been the recommended substrate for IFA. Monkey esophagus is used for the detection of both basement membrane antibodies and intercellular substance antibodies.<sup>1,5</sup> The intercellular substance antibody has been associated with the presence of a variety of disorders of the skin.<sup>3</sup> The detection of the basement membrane antibody has been associated with the presence of a variety of bullous pemphigoid autoimmune disorders of the skin.<sup>3</sup>

### Principles:

The primary test reaction involves circulating anti-epidermal antibodies present in the patient's serum, which attach to their homologous epidermal antigens. This occurs during the incubation period while the serum covers the antigen surface. A secondary reaction then follows a rinsing period which removes all unbound human antibody. The reagent used in the secondary reaction is a fluorescein labeled anti-human globulin conjugate containing Evans Blue Counterstain. The antigen surface is then thoroughly rinsed free of unbound conjugate and viewed under the appropriate fluorescent microscope for various morphological patterns of epidermal fluorescence which can be visually identified.

### Materials Provided:

#### Storage & Stability of Components:

- FITC IgG Conjugate No. AD CGEM2 (3.0 ml) (for use with Primate Substrates) is to be stored at 2-8 C upon receipt. The conjugate is stable at this temperature until expiration date on the vial label.
  - The antigen slides of monkey esophagus sections must be stored at 2-8 C upon receipt. Check label for specific expiration date.
  - ASA positive control No. AD PCBM (1.0 ml) for Basement Membrane reaction should be stored at 2-8 C upon receipt. Check label for specific expiration date.
  - ASA positive control No. AD PCIC (1.0 ml) for Intercellular Substance reaction should be stored at 2-8 C upon receipt. Check label for specific expiration date.
  - Universal negative control No. AD NC (1.0 ml) should be stored at 2-8 C upon receipt. Check label for specific expiration date.
  - Buffer Pack No. AD PBS1 - Phosphate Buffered Saline is stable at room temperature storage. Check label for specific expiration date. The reconstituted Buffer does not contain preservatives and should be stored at 2-8 C. Care should be taken to avoid contamination.
  - Mounting Medium No. AD TMM3 is stable when stored at 2-8 C. Check label for specific expiration date.
- Note: All kit components are available separately. Please see the current ALPHADIA Corporation Catalog for more details.

### Additional Materials Required but not Provided:

Test tubes and rack or microliter system  
Disposable pipettes  
Staining Dish and Slide Forceps  
Moisture Chamber  
Volumetric Flask (500 ml)  
Distilled H<sub>2</sub>O  
Fluorescence Microscope  
Paper Towels - lint free

### Reagent Preparation:

- Buffer Pack No. AD PBS1. Rehydrate buffer with 1 liter of sterile distilled water.

### Specimen Collection:

Serological specimens should be collected under aseptic conditions. Hemolysis is avoided through prompt separation of the serum from the clot. Serum should be stored at 2-8 C if it is to be analyzed within a few days. Serum may be held for 3 to 6 months by storage at -20 C or lower. Lipemic and strongly hemolytic serum should be avoided. When specimens are shipped at ambient temperatures, addition of a preservative such as 0.01% thimerosal or 0.095% sodium azide is strongly recommended.

### Test Instruction:

**Screening:** dilute test serums 1:20 in PBS. **Titration:** set up doubling dilutions of serum starting at 1:20, 1:40, 1:80, 1:160, etc.

- Once slides reach room temperature tear slide envelope at notch. Carefully remove the slide and avoid touching the antigen areas. The slide is now ready to use.
- Place a drop of diluted serum (20 to 30 µl) and controls over the antigen wells.
- Place slide with patient's serum and controls in a moist chamber for 30 minutes at room temperature (approximately 24 C).
- Remove slide from moisture chamber and allow the serum to run off onto a piece of paper towel. Using a wash bottle, gently rinse remaining sera from slide being careful not to aim the rinse stream directly on to the well.
- Wash in PBS for five minutes. Repeat using fresh PBS.
- Place a blotter on the lab table with absorbent side up. Remove slide from PBS and invert so that tissue side faces absorbent side of blotter. Line up wells to blotter holes. Place slide on top of blotter. **Do not allow tissue to dry.** Wipe back of slide with dry lint free paper towel. Apply sufficient pressure to slide while wiping to absorb buffer.
- Deliver 1 drop (25-30 µl) of conjugate per antigen well. Repeat steps 3-6.
- Place 4-5 drops of mounting medium on slide.
- Apply a 22 x 70 mm coverslip. Examine the slide under a fluorescent microscope. Note: To maintain fluorescence, store mounted slide in a moisture chamber placed in a dark refrigerator.

### Quality Control:

- Positive and negative serum controls must be included in each day's testing to confirm reproducibility, sensitivity and specificity of the test procedure.
- The negative serum control should result in little (+) or no fluorescence. If this control shows bright fluorescence, either the control, antigen, conjugate or technique may be at fault.
- The positive serum control should result in bright 3+ to 4+ fluorescence. If this control shows little or no fluorescence, either the control, antigen, conjugate or technique may be at fault.
- In addition to positive and negative serum controls, a PBS control should be run to establish that the conjugate is free from nonspecific staining of the antigen substrate. If the antigen shows bright fluorescence in the PBS control repeat using fresh conjugate. If the antigen still fluoresces, either the conjugate or antigen may be at fault.

### Results:

- Diffuse staining throughout the tissue is considered non-specific and should be considered a negative result.
- Staining of the basement membrane (BM) of the epidermis is considered positive and is associated with 70% of bullous pemphigus cases.
- Staining of the intercellular substance (ICS) of the prickle cell layer of the epidermis is considered positive 90% of pemphigus cases.
- The titer is the highest dilution of the patient's serum, showing a weak 1+ fluorescence of the ICS or BM. Titers of 1:20 or greater are clinically relevant for both patterns.

### CORRELATION OF BULLOUS DISEASE TO DIAGNOSTIC ASSAYS

Diseases	Biopsy Findings Direct	Serum Findings Indirect	Relevance of Immunofluorescence Findings
	Immunofluorescence	Immunofluorescence	
Pemphigus, all forms	ICS deposits of IgG	ICS antibodies	DIF and/or IIF diagnostic
Pemphigus, all forms	BM linear IgG and/or C also sometimes other Ig and F	BM antibodies IgG about 70%	DIF or IIF diagnostic
Scarring pemphigoid and Brunsting-Perry variety	BM linear IgG and/or C; also other Ig-90%	BM antibodies IgG about 20%	DIF or IIF diagnostic
Herpes gestations	BM linear, mostly c; also Ig	BM antibodies IgG about 20%	DIF or IIF diagnostic
Heratitis Herpetiformis typical	IgA also F and C granular or fibrillar in dermal papillae 85%		DIF diagnostic
Linear IgA dermatosis	IgA linear; 15% at BM; also F & C	No specific antibodies or BM antibodies IgA	DIF diagnostic

### KEY TO ABBREVIATIONS:

ICS = Inter-Cellular Substance  
BM = Basement Membrane  
C = Complement  
DIF = Direct Immunofluorescence  
IIF = Indirect Immunofluorescence

### Limitations of Procedure:

- No diagnosis should be based on a single serologic test since various host factors must be taken into consideration.
- Patients with Lydel's toxic neurolysis, extensive burns and myasthenia gravis may demonstrate intercellular substance staining.
- Additional confirming tests for Bullous diseases are skin biopsy, for direct immunofluorescent analysis and electron microscopy study.

### Precautions:

- All human components have been tested by radioimmunoassay for (HB<sub>s</sub>A<sub>g</sub>) and HTLVIII/LAV by an FDA approved method and found to be negative. (Not repeatedly reactive). However, this does not assure the absence of HB<sub>s</sub>A<sub>g</sub> or HTLVIII/LAV. All human components should be handled with appropriate care.
- The sodium azide (0.095%) included in the controls and conjugate is toxic if ingested.
- Do not use components beyond their expiration date.
- Follow the procedural instructions exactly as they appear in this insert to insure valid results.
- For *in vitro* diagnostic use.
- Handle slides by the edges since direct pressure on the antigen wells may damage the antigen.
- Once the procedure has started do not allow the antigen in the wells to dry out. This may result in false negative test results, or unnecessary artifacts.
- Use separate pipette tips for each sample and reagent to avoid cross contamination.
- Reagents should be inspected for evidence of bacterial or fungal contamination.
- Do not reuse substrate slide.


Component	AD PBS1 PBS Powder Packets AD TMM3 Mounting Medium	<b>Precautionary Statement Prevention:</b> P264 Wash thoroughly after handling. P280 Wear protective gloves and clothing.
Pictogram		<b>Response:</b> P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if easy to do. Continue rinsing.
Signal Word	<b>WARNING</b>	
Hazard Statement	H319 Causes serious eye irritation	P337+P313 If irritation persists, get medical advice/attention.

**Bibliography:**

1. Beutner et al.: Labeled antibody studies. In: Immunopathology of the Skin. Dowden, Hutchinson and Ross. 1973.
2. Fry and Beah: Immunological Aspects of Skin Diseases. John Wiley and Sons. 1974.
3. Beutner et al.: Autosensitization in Pemphigus and Biliious Pemphigoid. Charles C. Thomas. 1970.
4. Nakamura and Deodhar: Laboratory Tests in the Diagnosis of Autoimmune Disorders. Am Soc Clin Pract, Chicago, Illinois, p. 103-11 1976.
5. Beutner et al.: Immunopathology of the Skin, Second Edition, John Wiley and Sons. 1976.

**ALPHADIA sa/nv**  
DIAGNOSTIC PRODUCTS  
Avenue Vésale 26  
B1300 WAVRE  
BELGIUM  
TEL : 32 (0) 10 24 26 49  
FAX : 32 (0) 10 24 55 99  
contact@alphadia.be

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