## 1. INTENDED USE

The VIR-ELISA CMV-IgG AVIDITY is a supplementary reagent kit for the determination of CMV-specific IgG avidity in human serum and plasma. It has to be used in combination with the VIR-ELISA ANTI-CMV-IgG [REF] EG 101.

[REF] EAv 101 ▽ 48 [IVD]

## 2. PRINCIPLE OF THE ASSAY

The detection of antibodies is based on the principle of an enzyme-linked immunosorbent assay (ELISA). The purified, homogeneous antigen is fixed to each well of the microtiterstrips. Any specific antibodies present in the patient's sample (set up in duplicate) are bound during the first incubation.

After removing unbound material by washing, avidity reagent [UREA] is added to one well of the duplicate and wash solution (ready to use) to the other well of the duplicate. Using the special reagent for avidity the binding capacity of immune complex becomes dissolved but only for antibodies with low affinity whereas antibodies with high affinity keep associated with the antigen.

After removing unbound material by a further wash cycle, the presence of specific antibodies which are still bound to the antigen are detected using Anti-human IgG during the third incubation.

Excess peroxidase conjugate is then removed and TMB substrate is added, resulting in the development of a blue colour. The enzyme reaction is terminated by the addition of a stop solution. The intensity of the yellow colour thus developed is proportional to the concentration of antibodies in the sample.

## 3. DIAGNOSTIC RELEVANCE

Human cytomegalovirus (CMV) is a herpesvirus which is ubiquitously distributed in the human population. CMV is the most common cause of congenital infection, occurring in approximately 1% of all live births. Since CMV infections in immunocompetent individuals and pregnant women are asymptomatic or accompanied by symptoms not specific for CMV, laboratory methods are needed to diagnose CMV infection. In the absence of seroconversion, CMV-specific immunoglobulin M (IgM) is a sensitive and specific indicator of active or recent CMV infection.

However, the presence of CMV IgM is not a specific indicator of primary CMV infection as it is often produced during nonprimary infections. The measurement of the CMV IgG avidity index has been shown to be useful in identifying and excluding primary CMV infections with no pregestational CMV serology. Avidity describes the binding strength of specific antibody to antigen. It was found to be low in the first phase after primary infection but then to increase over time.

Detection of low-avidity CMV IgG in specimens indicates that primary CMV infection has occurred within the past 6 to 8 months, whereas detection of high-avidity CMV IgG excludes primary infection.

## 4. COMPONENTS OF THE TEST KIT

The test kit includes:

### AVIDITY REAGENT [UREA]

One vial containing 5,0 ml of ready-to-use avidity reagent. The avidity reagent contains 0,095% sodium azide as preservative.

### POSITIVE CONTROL / LOW AVIDITY [CONTROL LA]

One vial of 1.2 ml, containing human serum (low avidity) with 0,095% sodium azide as preservative. Ready to use.

*The safety data sheet (MSDS) is available upon request.*

## 5. STORAGE AND STABILITY

Store all reagents at 2-8°C. Protect them from intense light and do not freeze. The expiration date of each component is indicated on the respective vial label. Do not use reagents beyond the expiration date.

## 6. MATERIALS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

- VIR-ELISA ANTI-CMV-IgG [REF] EG 101
- Test tubes for sample dilution
- Timer
- Micropipettes, multipipettes 10-1000 µl
- One-liter graduated cylinder, distilled water
- ELISA washer or multichannel pipette
- Spectrophotometer for micro-plates (450 nm/ reference wavelength 630/620 nm)
- Paper towels, pipette tips

## 7. WARNINGS OR PRECAUTIONS SAFETY PRECAUTIONS

The ELISA test is for [IVD] use only.

1. Only qualified and well-trained employees should carry out the assay procedure.
2. The instruction for use describes the applicable test method. In case of modification or applications others than the intended use, or the use of automatic processors, the user has to validate the procedure and take the responsibility for it.
3. Seal all bottles properly after use in order to avoid bacterial contamination. All samples and kit components should be considered potentially infectious. All control samples have been tested for Hepatitis Bs antigen, anti-HIV I and II, anti-HCV (CE/FDA) and found to be negative.
4. Avoid contact with skin and mucous membranes when handling reagents, which contain preservatives (see kit contents). Wash thoroughly with water in case of contact and possibly look up a doctor.
5. Reagents containing sodium azide may react with lead and copper plumbing, building up explosive metal acids. Flush with sufficient water when disposing of reagents.
6. For disposal the legal regulations have to be followed.

## 8. SPECIMEN COLLECTION AND STORAGE

1. Microbially contaminated specimen may cause interference.
2. Lipaemic, hemolytic or icteric samples should only be tested with reservations although in our testing no negative influence has been found.
3. Suitable specimens are serum or plasma (heparinized, EDTA) samples obtained by standard laboratory techniques.
4. The samples should not be heat-inactivated since non-specific results may occur.
5. Patient samples should be stored at  $\pm$  2-8°C. For long term storage  $\pm$  -20°C or lower is recommended. Avoid repeated freeze-thaw cycles.
6. **Note:** Diluted patient samples must be used on the same day.

## 9. GENERAL INFORMATION **Please note: important !**

### Impact of Antibody Concentration on the Determination of Avidity

Due to the competition of specific antibodies with low and high affinity to the same binding site of corresponding antigen (solid phase on microtiter plate), the concentration of antigen-specific antibodies will influence the level of the avidity index. In order to keep the ratio between antigen and antibody constant, the antibody concentration of the sample has to be within the linear measuring range, this means it has to be > than the upper value of the cut-off range and < than the value of the positive control ([CONTROL+]) from test kit [REF] EG 101. The avidity indices of samples with an antibody concentration within this range can be calculated right away without further dilution of the samples. But samples with an antibody concentration less than 1,5 ratio might cause a false classification of the avidity index and should therefore be retested with a new serum sample. For the correct determination of avidity, samples higher than the positive control ([CONTROL+]) should be prediluted prior to the test, depending on the previously measured CMV-IgG concentration. The following dilutions are recommended:

Samples > than the [CONTROL+] and with a ratio >5,5 should be prediluted 1:10 and

Samples > than the [CONTROL+] and with a ratio < 5,5 should be prediluted 1:2.

If the prediluted samples lie again above the [CONTROL+], they must be diluted further (e.g. 1:20) and if they lie below the cut-off value, the samples must be diluted less (e.g. 1:5).

## 10. REAGENT PREPARATION

The determination of CMV-specific IgG avidity is carried out under utilisation of the VIR-ELISA ANTI-CMV-IgG [REF] EG 101.

**Bring all reagents to room temperature (21°C-25°C) prior to use !**

[WASHBUF] : Dilute the [WASHBUF25x] 1:25 with distilled water e.g. add 40ml of [WASHBUF25x] to 960 ml distilled water and mix well.

**Dilution of samples IgG:** Dilute patient samples 1:101 with [SPE DIL] e.g. 10µl sample + 1ml [SPE DIL] or resp. 10µl prediluted sample (e.g. 1:2) with 1ml [SPE DIL]; mix thoroughly.

**Controls are ready to use.**

Take the required [MTS] out of the foil packets and place them in the holder. Possibly remaining wells of a [MTS] have to be stored at 2-8°C tightly sealed in the plastic bag provided, with the desiccant inside.

## 11. PIPETTING AND INCUBATION STEPS

- A. Pipette 100 µl of sample diluent into well A1 (Blank). Pipette 100µl of the controls or diluted patient sample into the wells. The **CONTROL+**, the **CONTROL** and the patient samples must be set up in duplicates.
- B. Incubate the wells at room temperature (21...25°C) for 30 minutes, protected from intense light.
- C. Wash the wells four times as described in section N. **WASHING PROCEDURE**
- D. Add 100µl of ready to use avidity reagent **UREA** to one well of the duplicate, add 100µl of the diluted wash solution to the other well of the duplicate and to all other wells (Blank, negative control and cut-off-control).
- E. Incubate the wells at room temperature (21...25°C) for exactly 5 minutes, protected from intense light.
- F. Wash the wells four times as described in section k. **WASHING PROCEDURE**
- G. Add 100µl of ready-to-use peroxidase conjugate to each well.
- H. Incubate the wells at room temperature (21...25°C) for 30 minutes, protected from intense light.
- I. Repeat washing as in section C above.
- J. Add 100µl of ready-to-use TMB substrate to each well.
- K. Incubate the wells at room temperature (21...25°C), in the dark for 10 minutes
- L. Add 100µl of stop solution to each well. Tap gently to ensure homogenous color distribution and read within 10 minutes.
- M. To read the plate, make sure that the bottom is free from moisture and that no air bubbles are in the wells. Read the absorbance of the well contents at 450nm on a suitable plate reader. On readers equipped with a dual wavelength facility set the reference filter to 620/630 nm.

### ATTENTION:

The absorbance (OD) of the Blank must be always subtracted from the OD values of the controls and samples.

### PROCEDURAL NOTES

Do not allow the wells to dry out between incubations.  
Comply with the given incubation temperatures and times.

### N. **WASHING PROCEDURE**

The washing procedure can be done manually with a multichannel pipette or on an automatic plate washer. Empty the wells, invert and tap dry on paper towel. Wash four times with a soaking time of approx. 30 seconds (300 µl ).

## 12. CALCULATION OF RESULTS

For each sample and positive control, calculate the percent ratio between the O.D. of the well treated with Avidity Reagent and the O.D. of the well treated with Wash Solution and multiply with the factor 100. This ratio will be the sample/control percent Avidity.

$$\frac{\text{O.D. value (with avidity reagent)}}{\text{O.D. value (with wash solution)}} \times 100 = \text{Avidity Index (AI \%)}$$

Interpretation of sample results:

RESULT	DEFINITION
low avidity IgG	avidity index < 40%
equivocal range	avidity index 40% - 45%
high avidity IgG	avidity index > 45%

Determination of CMV-specific IgG antibody avidity could be used as a diagnostic marker to distinguish between primary infection and non-primary infection. An avidity index of less than 40% indicates a primary infection acquired within the past 6 – 8 months. The range between 40% and 45% is defined as equivocal. Equivocal results should be retested. If the equivocal result is confirmed it is not possible to determine the time of infection. The presence of high avidity IgG, that means an avidity index >45%, essentially excludes the possibility that infection occurred within the previous 4 months.

## 13. VALIDITY OF THE ASSAY

All controls should be carried out with every test run.

The test must comply with the following validation criteria:

- Avidity index of the **CONTROL** should be <25%.
- Avidity index of the **CONTROL+** should be >50%.

See also validation criteria of the VIR-ELISA ANTI-CMV-IgG **REF** EG 101.

If controls give invalid levels then results from the test samples are invalid too and retesting is required.

## 14. PERFORMANCE CHARACTERISTICS

### Diagnostic validity

#### Clinical specificity

The clinical specificity of the test was evaluated testing 166 anti-CMV-IgG positive samples from blood donors, pregnant women and patients. All the samples showed an avidity index >45% (mean value = 65%). **Specificity value = 100 %**. In addition the clinical specificity of the test was evaluated testing 18 samples coming from 9 patient panels with a history of CMV reinfections/ CMV past infections. All the samples showed an avidity index >45% (mean value = 63%). **Specificity value = 100 %**.

#### Clinical sensitivity

The clinical sensitivity of the test was evaluated testing two commercially available seroconversion panels. These panels, one a 15 member seroconversion panel and the other a 25 member seroconversion panel, depict the onset and decline of IgM and IgG antibodies over a period of time. The results of the avidity testing have been compared to the results of a reference assay (ELISA). All the samples showed an avidity index <40% (mean value = 21%) The sensitivity was calculated based on the results found. **Sensitivity value = 100 %**. In addition the clinical sensitivity of the test was evaluated testing 24 samples coming from 16 patient panels with evidence of acute primary CMV infections. All the samples showed an avidity index <40% (mean value = 23%). **Sensitivity value = 100 %**.

The results refer only to the groups of samples investigated.

### Precision and reproducibility

Intra-assay reproducibility was determined with samples of different avidity by testing each sample for at least 22 times in one test run. The calculated coefficients of variation (CV) of the samples were < 10%.


Inter-assay reproducibility was determined by testing samples of different avidity in 10 independent test runs. The calculated CV's of the samples were < 10%.

## 15. **Bib**

1. Avidity of Immunoglobulin G Directed against Human Cytomegalovirus during Primary and Secondary Infections in Immunocompetent and Immunocompromised Subjects. Clinical and Diagnostic Laboratory Immunology, July 1997, p. 469-473, Vol. 4, No. 4.
2. A New Method with General Diagnostic Utility for the Calculation of Immunoglobulin G Avidity. Clinical and Diagnostic Laboratory Immunology, Sep. 1999, p. 725-728, Vol. 6, No. 5.
3. Diagnosis of and Screening for Cytomegalovirus Infection in Pregnant Women. Journal of Clinical Microbiology, Sep. 2005, p 4713-4718, Vol. 43, No. 9.

**Viro-Immun Labor-Diagnostika GmbH**, In der Au 29, D-61440 Oberursel, Germany -Tel.: +49-6171-6281-00 - Fax: +49-6171-6281-12  
email: info@viro-immun.de

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