# VIR-ELISA ANTI - C. PNEUMONIAE- IgM

## 1. INTENDED USE

Enzyme immunoassay for qualitative determination of IgM specific antibodies to Chlamydia pneumoniae in human serum and plasma

IgM [REF] EM 128 ▼ 96 [IVD]

## 2. PRINCIPLE OF THE ASSAY

The detection of antibodies is based on the principle of an enzyme-linked immunosorbent assay (ELISA). The purified, homogeneous antigen is fixed to each well of the microtiterstrips. Any specific antibodies present in the patient's sample are bound during the first incubation. After removing unbound material by washing, the presence of specific antibodies is detected using Anti-human IgM conjugate during the second incubation.

Excess peroxidase conjugate is then removed and TMB substrate is added, resulting in the development of a blue colour. The enzyme reaction is terminated by the addition of a stop solution. The intensity of the yellow colour thus developed is proportional to the concentration of antibodies in the sample.

## 3. DIAGNOSTIC RELEVANCE AND INTERPRETATION OF RESULTS

The diagnosis of *C. pneumoniae* infections encloses the isolation of the pathogen from cell culture, direct antigen detection, nucleic acid amplification tests, and serology. Difficulties in sample collection and inaccessibility of culture and non-culture methods have made serology the method of choice. Besides the microimmune fluorescence method (MIF), antibodies to *C. pneumoniae* can also be detected by Enzymimmunoassay.

Primary chlamydial infection is characterized by a predominant IgM response within 2 to 4 weeks and a delayed IgG and IgA response within 6 to 8 weeks. After acute *C. pneumoniae* infection, IgM antibodies are usually no longer detectable within 2 to 6 months. IgG antibody titres usually decrease slowly, whereas IgA antibodies tend to disappear rapidly. When primary chlamydia infection is suspected, the detection of IgM is highly diagnostic. However, in recurrent or chronic infections the prevalence of IgM is low and therefore absence of IgM does not necessarily exclude on-going infection. In reinfections, IgG and IgA levels rise quickly. IgA antibodies have shown to be a reliable immunological marker of primary, chronic and recurrent infections. These antibodies usually decline rapidly to baseline levels following treatment and eradication of the chlamydia infection. The persistence of elevated IgA antibody titres is generally considered as a sign of chronic infection.

IgG	IgM	IgA	Interpretation	Recommendation
-	-	-	No specific antibodies detectable	No indication of infection. By suspect of acute infection control tests within 10 to 14 days are recommended. In case of justified clinical suspicion direct antigen detection should be performed.
-	+	+	Possibility of an early stage of infection or solitary, persisting IgA	Retest IgG, IgM and IgA after 10 to 14 days to determine seroconversion.
+	+	+	Acute or recent infection possible;	<b>Monitoring</b> of IgG, IgM and IgA antibodies (Sample collection within 10-14 days):
+	-	+	Reinfection or chronic infection possible	
-	+	-	Recent infection possible	<b>Monitoring</b> of IgG, IgM and IgA antibodies to determine seroconversion.
+	-	-	Indication of a past infection; Reinfection or chronic infection might also be possible	In case of clinical suspicion control tests of IgG and IgA antibodies within 10 to 14 days are recommended.

- negative, + positive

### Note

In cases of fresh acute chlamydial infections the serological antibody results may be negative despite clinics and positive antigen detection. If a serological confirmation of a positive antigen result is desired we recommend testing after 10-14 days for seroconversion.

## 4. PERFORMANCE CHARACTERISTICS

### Specificity / Sensitivity

135 samples were tested parallel in VIR-ELISA ANTI-C. PNEUMONIAE-IgM and comparison methods (ELISA/MIF). The specificity and sensitivity are based on the results found.

**Specificity:** IgM 96,8% **Sensitivity:** IgM 81,8%

For the calculation of the **specificity and sensitivity** of the ELISA test, equivocal results were defined as positive results. The results refer to the groups of samples investigated.

### Precision and reproducibility

Intra-assay reproducibility was determined by testing samples of different levels of antibody reactivity for at least 22 times in one test run. The coefficient of variation (CV) of the reactive IgM samples was < 10%.

Inter-assay reproducibility was determined by testing samples of different levels of antibody reactivity in 10 different test runs. The CV of the reactive IgM samples was < 10%.

### Cross-reactivity

Crossreactions between *C. pneumoniae* and other members of the chlamydiae family (*C. trachomatis* and *C. psittaci*) may occur due to the similarity of the antigens of these bacteria. In the case of positive or equivocal results an appropriate test (MIF) is advisable to exclude the above-mentioned species.

In addition, cross reactivity with antinuclear antibodies and heterophilic antibodies cannot be excluded in individual cases.

## 5. [Bib]

- Diagnostik von Chlamydien-Infektionen, J. Lab Med 2004; 28 (2):144-153
- Epidemiologisches Bulletin des Robert Koch Institutes; 6. April 2001 / Nr. 14 „Erkrankungen durch Chlamydia psittaci und Chlamydia pneumoniae: Ratgeber Infektionskrankheiten 23. Folge“.
- Serodiagnostics of Chlamydial Infections-Significance of Positivity in IgA and/or IgM Antibody Classes only. P. Hejnar, D. Koukalová; Biomed. Papers 146(2), 33- 35 (2002)

## 6. KIT COMPONENTS

- MICROTITERSTRIPS** [MTS]  
One microtiterplate is supplied which contains 12 microtiterstrips of 8 breakpart wells. The wells are coated with purified, inactive antigen. Strips are colour-coded.
- PEROXIDASE CONJUGATE** [CONJ|POD|IgM]  
One vial containing 12 ml of ready-to-use sheep anti-human IgM Peroxidase conjugate. Peroxidase conjugate contains 0.049% Thimerosal as preservative.
- NEGATIVE CONTROL** [CONTROL-]  
One vial of 1.2 ml containing human serum with 0,095% sodium azide as preservative. Ready to use.
- CUT OFF CONTROL** [CUTOFF]  
One vial of 1.2 ml containing human serum with 0,095% sodium azide as preservative. Ready to use.
- POSITIVE CONTROL** [CONTROL+]\*  
One vial of 1.2 ml containing human serum with 0,095% sodium azide as preservative. Ready to use.
- IgG-RF-SORBENT** [RF-SORB]  
One vial of 7,5 ml containing precipitating Anti-Human IgG reagent for the elimination of IgG and rheumatoid factors from human sera with 0,095% sodium azide as preservative. Ready to use.
- TMB SUBSTRATE** [SUBS|TMB]  
One vial containing 13 ml of ready-to-use tetra-methylbenzidine (TMB) substrate.
- SAMPLE DILUENT** [SPE|DIL]  
One bottle containing 100 ml (or 2x50 ml) of ready-to-use sample diluent buffer. The buffer includes 0.05% microbial preservative.
- WASH SOLUTION 25X** [WASH|BUF|25x]  
One bottle containing 80 ml (or 2x40 ml) of wash solution concentrate.
- STOP SOLUTION** [SOLN|STOP]  
One bottle containing 15 ml of 0.95N H<sub>2</sub>SO<sub>4</sub> stop solution. Ready to use.
- ADHESIVE SEALS**  
**The safety data sheet (MSDS) is available upon request.**

\* Is no titre indicated on the label, the [CONTROL+] is a positive reference sample and has a ratio ≥ 1,500

## 7. STORAGE AND STABILITY

Store all reagents at 2-8°C. Protect them from intense light and do not freeze. The expiration date of each component is indicated on the respective vial label. Do not use reagents beyond the expiration date.

After opening, [MTS] must be stored at 2-8°C in the plastic bag with desiccant and are stable up to 4 weeks.

The diluted [WASHBUF] is stable up to 4 weeks if stored at 2-8°C.

Use only [MTS] with an intact vacuum packaging.

## 8.

### MATERIALS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

- Test tubes for sample dilution
- Timer
- Micropipettes, multipipettes 10-1000 µl
- One-litre graduated cylinder, distilled water
- ELISA washer or multichannel pipette
- Spectrophotometer for micro-plates (450-nm/ reference wavelength 630/620 nm)
- Paper towels, pipette tips
- Incubator for 37°C with high relative humidity

## 9. WARNINGS OR PRECAUTIONS SAFETY PRECAUTIONS

The ELISA test is for **IVD** use only.

1. Do not mix lot specific reagents, such as **MTS**, controls and **CONJ/POD** from different kit lots. The **SUBS/TMB** must not be from the original test kit, but the lot of the **SUBS/TMB** must be the same as indicated on the kit label. The **SPE/DIL** (except Immunocapture assays), the **WASHBUF/25x** and the **SOLN/STOP** can be used for all ELISA tests.
2. Seal all bottles properly after use in order to avoid bacterial contamination. All samples and kit components should be considered potentially infectious. All controls have been tested for Hepatitis Bs antigen, anti-HIV I, anti-HIV II, anti-HCV (CE/FDA), and found to be negative.
3. The **MTS** are coated with inactive antigen. However, normal laboratory precautions should be maintained when handling with infectious material. Do not pipette by mouth.
4. Avoid contact with skin and mucous membranes when handling reagents, which contain preservatives (see kit contents). Wash thoroughly with water in case of contact and possibly look up a doctor.
5. Controls containing sodium azide may react with lead and copper plumbing, building up explosive metal acids. Flush with sufficient water when disposing of reagents.
6. The **SOLN/STOP** (0.95 N H<sub>2</sub>SO<sub>4</sub>) contains sodium hydroxide which may irritate skin and mucous membranes. Wash thoroughly with water in case of contact.
7. For disposal the legal regulations have to be followed.

## 10. SPECIMEN COLLECTION AND STORAGE

1. Only qualified and well-trained employees should carry out the assay procedure.
2. The instruction for use describes the applicable test method. In case of modification or applications others than the intended use, or the use of automatic processors, the user has to validate the procedure and take the responsibility for it.
3. Microbial contaminated specimen may cause interferences.
4. Lipaemic, hemolytic or icteric samples should only be tested with reservations although in our testing no negative influence has been found.
5. Suitable specimens are serum or plasma (heparinized) samples obtained by standard laboratory techniques.
6. The samples should not be heat-inactivated since non-specific results may occur.
7. Results on tests using CSF are not available.
8. Patient samples should be stored at 2-8°C. For long-term storage -20°C or lower is recommended. Avoid repeated freeze-thaw cycles.
9. **Note:** Diluted patient samples must be used on the same day. In order to avoid interferences with high IgG titres and rheumatoid factors (RF) an absorption of IgG/RF with IgG-RF-Sorbent must be performed.

## 11. ASSAY PROCEDURE REAGENT PREPARATION

Bring all reagents to room temperature prior to use!

**WASHBUF:** Dilute the **WASHBUF/25x** 1:25 with distilled water  
e.g. add 40ml of **WASHBUF/25x** to 960 ml distilled water and mix well.

**Dilution of samples:** 10 µl sample + 500 µl **SPE/DIL** (sample dilution 1:51). Mix thoroughly!

**Absorption:** Mix diluted sample 1+1 with IgG-RF-Sorbent (sample dilution 1:101). Mix thoroughly!

e.g. for single determination: mix 75µl diluted sample (1:51) with 75 µl IgG-RF Sorbent.

Incubate for 30 minutes at room temperature (18...25°C) or incubate over night at 2° to 8°C.

For carrying out the test the corresponding amount (100 µl) can be taken directly out of the mixture without prior centrifugation.

**Controls are ready to use and they should not be pre-treated with IgG-RF-Sorbent.**

**Note:** To take into consideration pipetting time, it is recommended to repeat the **CUTOFF** after every 4 **MTS** (resp. after a pipetting time of >=5 min.) to evaluate the following patient's tests with the new calculated cut-off value.

Take the required **MTS** out of the foil packets and place them in the holder. Possibly remaining wells of a **MTS** have to be stored at 2-8°C tightly sealed in the plastic bag provided, with the desiccant inside.

## 12. PIPETTING AND INCUBATION STEPS

- A. Pipette 100µl of each control and each diluted patient sample into the wells. Pipette 100 µl of sample diluent into well A1 (Blank).
- B. Cover the **MTS** with an adhesive seal and incubate the wells for 60 minutes at 37°C (± 1°C) in an incubator with relative high humidity or in a humid chamber.
- C. Wash the wells four times as described in section k. **WASHING PROCEDURE**
- D. Add 100µl of ready-to-use peroxidase conjugate to each well.
- E. Cover the **MTS** with an adhesive seal and incubate the wells for 30 minutes at 37°C (± 1°C) in an incubator with relative high humidity or in a humid chamber.
- F. Repeat washing as in section C above.
- G. Add 100µl of ready-to-use TMB substrate to each well.
- H. Incubate the wells at room temperature (18...25°C), in the dark for 15 minutes
- I. Add 100µl of stop solution to each well. Tap gently to ensure homogenous colour distribution and read within 10 minutes.
- J. To read the plate, make sure that the bottom is free from moisture and that no air bubbles are in the wells. Read the absorbance of the well contents at 450nm on a suitable plate reader. On readers equipped with a dual wavelength facility set the reference filter to 620/630 nm.

**Attention:** The absorbance (OD) of the Blank must always subtracted from the OD values of the controls and samples.

### PROCEDURAL NOTES

Do not allow the wells to dry out between incubations.  
Comply with the given incubation temperatures and times.

### k. WASHING PROCEDURE

The washing procedure can be done manually with a multichannel pipette or on an automatic plate washer. Empty the wells, invert and tap dry on paper towel. Wash four times with a soaking time of approx. 30 seconds (300 µl).

## 13. SUGGESTIONS FOR TROUBLESHOOTING

In case that the ELISA instructions are followed strictly, the reagents are handled with care and the samples and reagents are pipetted carefully, the following kinds of errors can be avoided to a large extend.

ERROR	POSSIBLE CAUSES
No colorimetric reaction after addition TMB substrate	No Peroxidase conjugate dispensed, contamination of Peroxidase conjugate (possibly with control sera during pipetting) may cause an inactivation.
Generally too high reaction	Incorrect Peroxidase conjugate (i.e. not from original test kit), incubation time too long or incubation temperature too high, water quality for Washing Solution insufficient (low grade of deionisation)
Generally too weak reaction	Incorrect Peroxidase conjugate (i.e. not from original test kit), incubation time too short, incubation temperature too low.
Reagent blank too high	Incorrect pipetting of sample diluent, contaminated reagents, reagents expired, exceeding of incubation time and temperature, external contamination of the bottom of microtiterstrips, (clean carefully!).
False positive / negative samples	Incorrect dilution of samples, microbial contaminated specimen.
Unexplainable outliers	Contamination of pipettes, tips or containers or with metals (iron, copper etc.), insufficient washing.
High variation (within a series)	Reagents (including microtiterstrips) not pre-warmed to room temperature prior to use. Washer is not washing correctly!
High variation (from series to series)	Incubation conditions not constant (time, temperature) high variation of incubation temperature, controls and samples are not carried out at same time (same intervals) check pipetting order, person related variation, strips dried out after washing (results are not reproducible).

## 14. VALIDITY OF THE ASSAY

All controls should be carried out with every test run.

The test must comply with the following validation criteria:

- OD-value of the Negative Control should be  $< 0.100$ ,
- OD-value of the Cut-off Control should be  $> 0.200$ ,
- Ratio of the Positive Control/ Cut-off Value should be  $\geq 1.5$
- OD-value of the Blank should not be higher than  $0.100$ .

If controls give invalid levels then results from the test samples are invalid too and retesting is required.

## 15. CALCULATION OF RESULTS

### A QUALITATIVE CALCULATION

#### Calculation of "Cut-off Value"

The **Cut-off Value** is calculated from the absorbance of the Negative Control and the absorbance of the Cut-off Control and defines the Cut-off range.

$$\text{Cut-off Value} = \text{OD of the Negative Control} + \text{OD of Cut-off Control}$$

$$\text{CUT-OFF RANGE} = \text{CUT-OFF VALUE} \pm 10\%$$

Interpretation of sample results:

RESULT	DEFINITION
negative -	OD value sample $<$ Cut-off value <b>-10%</b>
equivocal	OD value sample $\geq$ Cut-off value <b>-10%</b> OD value sample $\leq$ Cut-off value <b>+10%</b>
positive +	OD value sample $>$ Cut-off value <b>+10%</b>

Equivocal results should be retested. Following the confirmation of the equivocal result the monitoring of the patient's antibodies is recommended in order to exclude unspesific reactions resp. cross-reactivity, which may also cause equivocal results.

### B CALCULATION OF RATIO (CUTOFF INDEX, COI):

Patient samples may also be quantified and interpreted by means of the calculation of the ratio (Cutoff Index, COI):

$$\text{COI} = \text{OD value of sample} / \text{Cut-off value},$$

whereby a ratio of 1.000 is equivalent to the Cut-off value.

Interpretation of sample results:






**Ratio  $> 1,1$**  positive result  
**Ratio  $< 0,9$**  negative result  
**Ratio  $0,9-1,1$**  equivocal result

Regarding diagnostic relevance and interpretation of results see **page 1**.

For further information please visit our website:

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