ZenTech

# ELIZEN Pneumococcus IgG Immunopotency level

ENGLISH (en)

FRANÇAIS (fr)





ITALIANO (it)

**ZenTech s.a.**  
Liège Science Park,  
Avenue du Pré-Ailly 10,  
4031 ANGLEUR, Belgium  
Tel. : + 32-(0)4-361.42.32  
Fax : + 32-(0)4-367.00.63

[info@zentech.be](mailto:info@zentech.be)

[www.zentech.be](http://www.zentech.be)

ISO15223	MEDICAL DEVICES SYMBOL	SYMBOLES APPLIQUES AUX DISPOSITIFS MEDICAUX	SIMBOLO DEI DISPOSITIVI MEDICALI
	Storage temperature limitation	Limites de températures	Limite temperatura di conservazione
	Lot number	Numéro de lot	Numero di lotto
	Expiration date	Date d'expiration	Data di scadenza
	Consult operating instructions	Lire les instructions	Consultare le istruzioni di funzionamento
	<i>In vitro</i> diagnostic device	<i>in vitro</i> diagnostic	Dispositivo per diagnostica <i>in vitro</i>
	Manufactured by	Fabriqué par	Fabbricato da
	Catalog number	Référence	Numero catalogo

	SYMBOLS (EDMA recommendations)	SYMBOLES (recommandations EDMA)	SIMBOLI (raccomandazioni EDMA)
	number of determinations (96)	Nombre de déterminations	Numero di determinazioni (96)
	Calibrators	Calibrateur	Calibratore
	Control serum	Sérum de contrôle	Siero di controllo
	Microtiterplate	Microplaque	Piastra per microtitolazione rivestita
	Enzyme conjugate	Conjugué HRP	Anti-hlgG-HRP coniugato
	Sample diluent	Diluant des échantillons	Diluyente per campioni
	Chromogen	Chromogène	Cromogeno
	Substrate buffer	Tampon Substrat	Buffer substrato
	Washing solution to be diluted ten-fold	Solution de lavage concentrée 10 fois	Soluzione di lavaggio da diluire 10 volte
	Blocking reagent	Solution d'arrêt	Reagente bloccante

# ENGLISH

## IMMUNO-ENZYME ASSAY FOR THE DETERMINATION OF ANTIBODIES ANTI-PNEUMOCOCCUS IN HUMAN SERUM AND IN PREPARATIONS OF IMMUNOGLOBULINS

E-DG-96

### CLINICAL APPLICATION

*Streptococcus pneumoniae* is a common host of the oro-pharyngeal flora. Diplococcus Gram positive, its capsule confers it its virulence and antigenicity. Antibodies against the polysaccharidic antigens of the capsule are protective of the infection. The rate of *Streptococcus pneumoniae* carriers, varies from age, season and environment infectious or not. Pneumococcus are the major cause of pneumonia acquired in communities and a frequent source of acute otitis of the middle ear, sinusitis and pharyngitis. Pneumonia, which are caused by pneumococcus, are reported to be the fifth cause of death in the United States (1). In Belgium, the incidence of worrying infections caused by pneumococcus is 20.000 cases/year (2). The rate of antibio-resistant strain grows constantly: in 1997, 10 % of strains resist to penicillin, 23 % to tetracyclin and 28 % to macrolides (3). The mortality rate is evaluated at 2.000 deaths/year.


If pneumococcus affect all individuals (70 % of healthy carriers), certain factors predispose to pneumococcus infection like age (below 5 years and over 60), broncho-pneumopathies, certain superior sero-digestive tract cancers, alcoholism, smoking and cardiac insufficiency. Hypo- or asplenyias represent a certain aggravating factor as well as hemopathies (lymphoma, multiples myeloma, child under 5 years drepanocytosis).

The assay allows the detection of antibodies anti-polysaccharidic capsule of the 23 most frequent serotypes (95 % isolated strains in Belgium) (3).

This method allows the detection of antibodies in serums containing low levels (concentration 1.000 times lower than the mean encountered by healthy individuals) and in hyper-immune serums obtained after vaccination.

### PRINCIPLE OF THE ASSAY

The method is based on immuno-enzymatic reaction. Patients serum prior diluted is added into the wells coated with polysaccharid capsule of 23 pneumococcus strains. After carefully washing of the microplate wells, horse radish peroxydase coupled to monoclonal anti-human IgG is added in each well. The monoclonal antibody anti-hlgG-HRP reacts with the complexes adsorbed on the solid phase. After washing, the chromogen solution is added in the wells and incubated. The reaction is stopped, the absorbance, proportional to the concentration of antibodies anti-pneumococcus, is measured at 405 and 450 nm with a spectrophotometer.

### MATERIEL PROVIDED

The ELIZEN Pneumococcus IgG kit contains the reagents needed to carry out 96 determinations.

All the reagents must be kept at 2-8°C. Expiration dates of each reagent are indicated on its labels.

#### 1. Coated microtiterplate

1 microtiterplate of 96 breakable wells coated with a mix of polysaccharid capsule highly purified of 23 types of pneumococcus. After opening, store in the minigrip? bag furnished with a dessicant.

#### 2. Calibrators

7 vials of lyophilized calibrators. These calibrators are produced from a pool of plasma of 6.000 donors. For the exact value, refer to the data sheet included.

Calibrators must be reconstituted with 1 ml of distilled water (stable 2 months after reconstitution). A trouble may appear after reconstitution. The trouble doesn't affect the performance of the assay.

Preservative: Merthiolate.

#### 3. Control Serum

1 vial of 1 ml lyophilized. Control serum must be reconstituted with 1 ml of distilled water (stable 2 months after reconstitution)

A trouble may appear after reconstitution. The trouble doesn't affect the performance of the assay.

Preservative: Merthiolate.

#### 4. Anti-hlgG-HRP conjugate

1 vial of 12 ml of mouse monoclonal antibody anti-human IgG coupled to horse radish peroxydase in Tris-HCl buffer containing bovine albumin, stabilizers and preservatives.

#### 5. Samples diluent

1 vial of 100 ml of Tris-HCl buffer containing casein, stabilizers and preservatives.

A trouble may appear in the buffer due to the presence of casein. This trouble doesn't affect the performance of the assay. The diluent should be thoroughly stirred before use.

Preservative: Merthiolate.

Ready for use.

Store at 2-8 °C. Stable until the expiration date written on the label.

#### 6. Washing solution

1 vial of 25 ml of phosphate buffer containing Tween 20.

If undissolved crystals are present, put them back in solution by placing the vial at 37 °C for a few minutes.

Bring to 250 ml with distilled water.

Preservative: Merthiolate.

The diluted wash solution may be kept at 2-8 °C for 1 month.

#### 7. Substrate buffer

1 vial of 15 ml of phosphate-citrate buffer containing H<sub>2</sub>O<sub>2</sub> (0,03%).

Before use, dilute with an equal volume of chromogen.

#### 8. Chromogen

1 vial of 15 ml of phosphate-citrate buffer containing DMSO and TMB.

Dilute with an equal volume of substrate buffer.

Avoid all exposure to direct light and use the prepared solution within one hour.

#### 9. Stop solution

1 vial of 15 ml of H<sub>2</sub>SO<sub>4</sub> 0,5 M.

### KIT REAGENTS

Reagents	Quantity	Physical state
Microtiterplate	1	Ready for use
Calibrators	7 x 1 ml	Lyophilized
Control serum	1 x 1 ml	Lyophilized
Anti-hlgG-HRP conjugate	1 x 12 ml	Ready for use
Samples diluent	1 x 100 ml	Ready for use
Wash solution	1 x 25 ml	Concentrated 10 x
Substrate buffer	1 x 15 ml	Ready for use
Chromogen	1 x 15 ml	Ready for use
Stop solution	1 x 15 ml	Ready for use

## MATERIALS REQUIRED BUT NOT SUPPLIED

- Adjustable automatic pipettes with disposable tips.
- Graduated cylinder.
- Distilled water.
- Aspirating pump or microplates washer.
- Spectrophotometer able to read at 450 ,405 and 620 nm.
- Dry Heater, adjustable at 37°C ± 1°C.

## WARNINGS AND PRECAUTIONS

For *in vitro* diagnostic use only.

**In order to obtain reproducible results, the following rules must be observed:**

- The assay must be performed after the instructions have been read and understood.
- Let the reagents reach room temperature before opening them.
- Do not mix reagents from different batches.
- Do not use reagents after their expiration dates.
- Use absolutely clean glassware.
- Use distilled water kept in absolutely clean tanks.
- Avoid contaminating the samples. To ensure this, use disposable tips for each sample and reagent.
- Residual amounts of sodium azide (NaN<sub>3</sub>) can destroy the conjugate's enzymatic activity.
- Traces of hypochlorite of soda can destroy many reagents' biological activity.
- Avoid contaminating the TMB and substrate buffer solutions with the conjugate or any other oxidant.
- Do not expose the reagents to direct light during storage or incubation.

**In order to avoid personal and environmental contamination, the following precautions must be observed :**

- Use disposable gloves while handling potentially infectious material and performing the assay.
- Do not pipette reagents by mouth.
- Do not smoke, eat, drink or apply cosmetics during the assay.
- Chromogenic substrate and Blocking Reagent should be handled with care. Avoid contact with skin, eyes and mucous membranes. In case of accident rinse thoroughly with running water.
- All material of human origin used for the preparation of this kit is tested negative for HBsAg, anti-HIV and anti-HCV. Since no test at present can guarantee complete absence of these viruses, all samples and reagents used for the assay must be considered potentially infectious; therefore, the assay waste must be decontaminated and disposed off, in accordance with established safety procedures. Disposable ignitable material must be incinerated; disposable non-ignitable material must be sterilized in autoclave for at least 1 hour at 121°C.
- Liquid wastes must be added with sodium hypochlorite at a final concentration of 3%. Let the hypochlorite act for at least 30 minutes. Liquid wastes containing acid must be neutralized with appropriate amounts of base before treating with sodium hypochlorite.
- Avoid splashing and aerosol formation; in case of spilling, wash carefully with a 3% sodium hypochlorite solution and dispose of this cleaning liquid as potentially infectious waste.
- Some reagents contain sodium azide as preservative; to prevent build-up of explosive metal azides in lead and copper plumbing, reagents should be discarded by flushing the drain with large amounts of water.

## SPECIMEN COLLECTION AND PREPARATION

It is recommended to use serum. Avoid using hemolyzed, lipemic or bacterially contaminated sera. Thoroughly mix thawed specimens before assay and avoid repeated freeze/thawing cycles which may cause loss of antibody activity and give erroneous results.

Samples may be stored at 2°- 8°C for up to 48 hours. For long term storage, samples should be stored frozen at -20°C or lower.

### Sample pre-dilution

Dilute samples 1/51 : 10 µl sample + 500 µl sample diluent  
Then realize from 1/51 dilution the following dilution :  
1/510 : 50 µl 1/51 dilution + 450 µl sample diluent

## ASSAY PROCEDURE

**Let the reagents stabilize at room temperature.**

- 1- Prepare the exact number of wells to perform the calibration curve and the samples dilutions in duplicate.
- 2- Add 100 µl of each calibrator , control serum and of each patients pre-diluted 1/510 in each well.
- 3- Incubate 1 hour at 37 °C.
- 4- Aspirate and wash 3 times with 300 µl of diluted washing.
- 5- Add 100 µl of conjugate anti-hlgG-HRP in each well.
- 6- Incubate 30 minutes at 37 °C.
- 7- Aspirate and wash 3 times with 300 µl of diluted washing.
- 8- Add 200 µl of the chromogen-substrate solution prepared extemporaneously in each well.
- 9- Incubate 15 minutes at 37 °C protected from light.
- 10-Add 100 µl of Blocking solution in each well at the same speed then the add of chromogen-substrate solution.
- 11-Read the absorbance of the wells (450, 405 and 620nm (as reference)). Reading must be completed within 20 minutes from the end of the assay.

## SCHEME OF THE ASSAY

Wells	Calibrator ( 0- 6)	Control serum	Diluted samples
<b>Reagents</b>			
<b>Calibrator ( 0- 6)</b>	100 µl	---	---
<b>Control serum</b>		100 µl	
<b>Diluted samples</b>	---	---	100 µl
- Incubate 1 hour at 37 °C.			
-Aspirate and wash 3 times with 300 µl of washing solution.			
<b>Conjugate anti-hlgG-HRP</b>	100 µl	100 µl	100 µl
-Incubate 30 minutes at 37 °C			
-Aspirate and wash 3 times with 300 µl of washing solution.			
<b>Solution substrate/chromogen</b>	200 µl	200 µl	200 µl
-Incubate 15 minutes at 37 °C and protecting from light.			
<b>Stop Solution</b>	100 µl	100 µl	100 µl
Read : 450 nm – 405 nm and 620 nm as reference.			

## CALCULATION AND INTERPRETATION OF THE RESULTS

In order to obtain a better sensitivity, the present method employs spectrophotometric reading at two wavelengths (450 and 405 nm).

For all O.D. overflow at 450nm, multiplied the O.D. 405nm by the correction factor calculated by the ratio between O.D. 450nm and O.D. 405nm. To establish this factor, used a sample reading at 405nm and 450nm.

Draw a calibration curve on millimetric graph paper, by plotting the calibrators concentration (x-axis) against the absorbance obtained for each calibrator (y-axis).


Corresponding anti- pneumococcus concentrations in mU/ml are obtained by interpolating the absorbances of each sample on the calibration curve.

## EXAMPLE OF CALCULATION

The values shown below must be considered as an example and must not be used in place of experimental data.

Description	O.D	B/Bmax (%)	Concentration (mU/ml)
Calibrator 0	0.043	0.8	0
Calibrator 1	0.282	5.5	125
Calibrator 2	0.639	12.4	250
Calibrator 3	1.329	25.8	500
Calibrator 4	2.693	52.3	1000
Calibrator 5	3.980	77.2	2000
Calibrator 6	5.154	100	4000
Control serum	3.165		1280
P1	0.061		25.1
P2	0.745		290
P3	2.735		1035
P4	4.702		2953

## TYPICAL CALIBRATION CURVE

## REFERENCE VALUES

It is recommended that each laboratory determines its own reference interval. Values reported below are only indicative.

The distribution of the values of anti-pneumococcus in a normal population is not gaussian and shows a wide dispersion (see the high standard deviation).

moyenne (mU/ml)	1181.6
écart-type (mU/ml)	1165.8
minimum (mU/ml)	3.0
maximum (mU/ml)	>4000
n	100

### Interpretation of results:

Lower than 250 mU/ml	Deficient in anti-pneumococcus antibodies
----------------------	---

## PERFORMANCES OF THE ASSAY

### SENSITIVITY

The sensitivity was calculated based upon the calibration curve and expressed as the minimal dose showing a significant difference from the Zero calibrator (mean value + 2 S.D.). This dose is 6.9 mU/ml.

### PRECISION

Precision was evaluated upon intra- and inter-assay variability, in 3 sera at different anti-pneumococcus concentrations.

#### Intra-assay

Sample	Mean (mU/ml)	S.D.	C.V. %	Replicates No.
a	260.5	22.9	8.8	20
b	1059.2	42.5	4.0	20
c	2600.1	193.2	7.4	20

#### Inter-assay

Sample	Mean (mU/ml)	S.D.	C.V. %	Assays No.
a	290.4	9.6	3.3	9
b	1026.2	27.7	2.7	9
c	2732.4	31.4	1.1	9

### ACCURACY

Accuracy of the method has been checked by the recovery and parallelism tests:

#### Recovery Test

One prediluted (1:510) sample, mixed with equal volumes of each calibration point, was tested.

Sample	Expected (mU/ml)	Measured (mU/ml)	Recovery %
S1	-	492	-
S1+ cal 0	246	286	116.2
S1+ cal 1	309	307	99.5
S1+ cal 2	371	343	92.5
S1+ cal 3	496	504	101.6
S1+ cal 4	746	741	99.4
S1+ cal 5	1246	1183	95.0
S1+ cal 6	2246	2379	105.9

#### Parallelism Test

One sera with high anti-pneumococcus concentrations was tested at different dilutions with the Zero calibrator.

Dilution	Expected (mU/ml)	Measured (mU/ml)	Recovery %
S1 undiluted	-	2691	100
1/2	1345	1472	109.4
1/4	673	681	101.2
1/8	336	324	96.4
1/16	168	114	68.0

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