# VIRO-BLOT ANTI – H. PYLORI - IgG/IgA

Instruction for use

## Immunoblot for qualitative determination of antibodies to *Helicobacter pylori* in human serum and plasma

IgG	<b>[REF]</b> BG 125/20	▽ 20	<b>[IVD]</b>	complete kit
IgG	<b>[REF]</b> BG 125/8	▽ 8	<b>[IVD]</b>	complete kit
IgA	<b>[REF]</b> BA 125/20	▽ 20	<b>[IVD]</b>	complete kit
IgA	<b>[REF]</b> BA 125/8	▽ 8	<b>[IVD]</b>	complete kit

### Antigen Characterisation

**Antigen:** *H. pylori* spezific antigen (acid glycine extract)

**Strain:** DSMZ 10242 (NCTC 11638)

**Purification:** Ultrasonication, centrifugation, filtration

**Cultivation:** Kelly-Medium

### Introduction

*Helicobacter pylori* ( previously *Campylobacter pylori* ) is a micro-aerophilic, Gram negative, curved motile pathogen. It was first cultured by Marshall and Warren in 1983 from patients with gastritis. *H. pylori* is adapted to survive in the gastric acid and to escape the local immune response. Bacterial urease of *H. pylori* creates a local alkaline micro-environment. In addition *H. pylori* may produce various toxins.

*H. pylori* is the etiologic agent of chronic antral type B gastritis and involved in the pathogenesis of gastroduodenal ulcer disease and gastric carcinoma. This bacterium is only found in close association with inflamed gastric epithelium, usually in the gastric antrum and the first part of the duodenum. About 95% of patients with duodenal ulcer and 70-80% of those with gastric ulcer have *H. pylori* infection. 50% of patients with Non-Ulcer-Dyspepsie are associated with *H. pylori* infection. Several studies showed an association between *H. pylori* colonization and the development of gastric carcinoma. In Western countries around 70% of gastric cancer and in developing countries around 80% of gastric cancer could be attributed to infection with *H. pylori*. *H. pylori* affects about 5 - 15% below the age of 20 years and about 65% in advanced years. Nowadays *H. pylori* is uncommon in young children. This theory is based upon the assumption of a decreasing rate of infection in younger age cohorts in recent decades due to a rise in living standards. In the developing countries the infection rate is higher due to worse hygienic situations and a low socio-economic status.

*H. pylori* is believed to be transmitted orally. Many researchers think that *H. pylori* is transmitted orally by means of fecal matter through the ingestion of waste tainted food or water. There is also evidence that spread by fecal oral contact with infected persons is likely. Therefore the prevalence of *H. pylori* infection within a family community is clearly more often.

### Principle of the WESTERN BLOT

The WESTERN BLOT may be used for the qualitative determination of human antibodies to *Helicobacter pylori* antigens.

For this Immunoblot an European strain of *Helicobacter pylori* (complete antigen) was treated with SDS and then electrophoretically separated in Polyacrylamide-gel according to molecular weight of the proteins. This was then transferred to Nitrocellulose (blotted) and cut into individual **[NTS]**. The **[NTS]** were then incubated with blocking buffer.

For the assay the **[NTS]** are incubated with the control and the samples to be investigated, antibodies present will bind specifically to the protein bands and form antigen-antibody complexes. Unbound serum or plasma is removed by means of washing. During the second incubation Anti-human IgG or IgA Conjugate binds to the antibody-antigen complex. In the final wash step excess Conjugate is removed and **[SUBS|CH]** substrate is added to the chamber.

The specific protein bands appear as purple bands. The reaction is then stopped by the addition of distilled water and the **[NTS]** may be dried.

With the aid of a kit specific evaluation diagram the band pattern may be assessed and evaluated.

The interpretation of the band pattern provides evidence of the antibody response.

### Materials required but not provided with the Kit

- Distilled water
- Pipettes
- Timer
- Measuring cylinder
- Tweezers
- Test tubes for sample- and conjugate dilution
- Water jet vacuum pump (to remove fluid from the incubation dish)
- Rocking platform or horizontal shaker

Kit contents ▽ 20 or ▽ 8

<b>[NTS]</b>	20 / 8	<b>Nitrocellulose test strips</b> with purified, separated <b><i>Helicobacter pylori</i></b> proteins
	3 / 1	<b>Incubation dish</b> with 8 chambers
<b>[WASHBUF 25x]</b>	25ml / 25ml	<b>Wash solution concentrate</b> (25x)
<b>[SPE DIL]</b>	2x21ml / 21 ml	<b>Dilution buffer</b> (with 0,05% microbial preservative) <b>[RTU]</b>
<b>[CONJ ALK 10x]</b> <b>[IgG   IgA]</b>	2,3 ml / 1,0ml	<b>Anti-human IgG or IgA Alkaline Phosphatase</b> (sheep) (10x) (preservative: 0,049% Thimerosal)
<b>[CONTROL - 10x]</b>	0,5 ml / 0,4ml	<b>Negative control</b> , human (preservative: 0,095% sodium azid) (10x)
<b>[SUBS CH]</b>	21ml / 9ml	<b>Chromogen-/ Substrat solution</b> <b>[RTU]</b>
	1	<b>Evaluation diagram with positive Control strip IgG</b>
	1	<b>Protocol sheet</b>

### Precautions

The WESTERN BLOT test is for **[IVD]** use only.

1. Do not mix lot specific reagents, such as **[NTS]**, Controls and **[CONJ|ALK|10x]** from different kit lots. The **[SUBS|CH]** doesn't have to be from the original test kit. But the lot of the **[SUBS|CH]** should be the same as indicated on the kit label. The **[SPE|DIL]**, the **[WASHBUF|25x]** can be used for all WESTERN BLOT tests irrespective of the lot and the nature of the WESTERN BLOT.
2. Seal all bottles properly after use in order to avoid bacterial contamination. All samples and kit components should be considered potentially infectious. All control samples have been tested for Hepatitis Bs antigen, anti-HIV I and II, anti-HCV and found to be negative.
3. The **[NTS]** are coated with inactive antigen. However, normal laboratory precautions should be maintained when handling with infectious material. Do not pipette by mouth.
4. Avoid contact with skin and mucous membranes when handling reagents, which contain preservatives (see kit contents). Wash thoroughly with water in case of contact and possibly look up a doctor.
5. Controls containing sodium azide may react with lead and copper plumbing, building explosive metal acids. Flush with sufficient water when disposing of reagents.
6. For disposal the legal regulations have to be followed.

### Storage and stability

Store all reagents at  $\pm$  2-8°C, protected from intense light. The expiration date of each component is indicated on the respective vial label. Do not use reagents beyond the expiration date.

The diluted **[WASHBUF]** is stable up to 4 weeks when stored at  $\pm$  2-8°C.

### Specimen collection and general handling

1. The assay procedure should be carried out only by qualified and well trained employees.
2. The instruction for use describes the applicable test method. In case of modification or applications others than the intended use, the user has to validate the procedure and take the responsibility for it.
3. Lipaemic, icteric, haemolysed or microbially contaminated specimen may cause interference.
4. Suitable specimens are serum or plasma (heparinized, EDTA, Citrat) samples obtained by standard laboratory techniques.
5. The samples should not be heat-inactivated since non-specific results may occur. Avoid repeated freeze-thaw cycles.
6. For ethical reasons, the use of CSF samples was only evaluated by testing pairs of CSF and serum samples (only if available), distributed by the "Institut für Standardisierung und Dokumentation" (INSTAND) in Germany. Therefore the use of CSF samples should be proved by testing a sufficient number of samples.
7. Patient samples should be stored at  $\pm$  2-8°C.  
For long term storage  $\pm$  -20°C or lower is recommended.
8. **Note:** Diluted patient samples must be used on the same day.

### Preparation of the test

Bring all reagents to room  $\pm$  \* prior to use !

\* room  $\pm$  up to maximally 25°C

**[WASHBUF]** : Dilute the **[WASHBUF|25x]** 1:25 with distilled water  
e.g. add 10ml of **[WASHBUF|25x]** to 240 ml distilled water and mix well.

**Control:** Dilute the Negative Control **1:10** with Specimen diluent.  
(900  $\mu$ l Specimen diluent + 100  $\mu$ l negative Control)

**Dilution of samples:** Dilute patient samples **1:101** with **[SPE|DIL]**  
e.g. 10 $\mu$ l sample + 1000 $\mu$ l **[SPE|DIL]**; mix thoroughly.

**Conjugate dilution:** Dilute the [CONJ|ALK|10x] during the first wash step with [SPE|DIL] e.g. for 5 strips: 4,5ml [SPE|DIL] + 0,5ml [CONJ|ALK|10x] = 5ml Conjugate dilution

**Please note:**

- Carefully remove the [NTS] using tweezers, holding the strip by the numbered end.
- Close the tube immediately after taking the [NTS].
- The [NTS] should be used in number order from low to high.
- Lay the required number of [NTS], one per chamber of the incubation dish, placing them with the number facing up.
- Ensure that the [NTS] remain covered with fluid and that they do not dry out between the incubation or wash steps.

The [CONTROL-] should be used for both IgA and IgG.

**WESTERN BLOT Assay procedure**

Incubation at room T \*

**First incubation (sample incubation) 60 min.**

Carefully pipette 1000 µl diluted Control or 1010 µl pre-diluted patient sample into the incubation chamber to coat the strip and incubate **60 min.** on the rocking platform.

**First wash step**

Completely remove patient samples (with the vacuum pump) and wash each strip **1 x** with **2ml** prepared [WASH|BUF] and immediately drain. Next wash the [NTS] a further **3 x** with **2 ml** [WASH|BUF] for **5 min.** each time. (Incubation on the rocking platform).

**Second incubation (conjugate incubation) 30 min.**

Completely remove the [WASH|BUF] and pipette **1 ml** pre-diluted [CONJ|ALK] in each chamber. Incubate for **30 min.** on the rocking platform.

**Second wash step**

(see first wash step)

**Third incubation (Substrate incubation) 15 min.**

Completely remove the [WASH|BUF] and pipette **1 ml** [RTU|SUBS|CH] in each well. Incubate from 10 min. to **15 min.** on the rocking platform. Remove the [SUBS|CH] completely from each compartment.

- Stop the enzyme-substrate reaction with **2x 2 ml** for 5 min. distilled water each.

Remove the [NTS] carefully from the compartment and dry completely on filter paper. The protein bands may be identified and evaluated with aid of the Evaluation diagram.

For the purposes of documentation the [NTS] should be kept dry and protected from the light.

**Test results**

Evaluation of the strips manuell:

Match the dry [NTS] at the edge of the protocol sheet and fix the [NTS] with tape or similar at the numbered ends of the [NTS].

The identification of the bands is possible with the positive control strip card. The immuno specific bands are signed on the positive control strip card.

The control card will be matched at the incubated fixed patient [NTS], the black line have to agree. Now it is possible to read directly the immuno specific bands.

Record the recognisable bands according to their molecular weight and intensity on the protocol sheet.

The control strip card is specific to each kit. All high specific and specific bands are signed on the control card.

Interpretation of result by using special Software ( Seroline, Viro-Immun)

The result of a [NTS] can be interpreted by using a special software. Therefore dried and stained [NTS] should be fixed onto the protocol sheet which is included in the manual. It is recommended to glue the [NTS] (near the number of the strip) at the indication line of the protocol sheet.

The identification of each immuno reactive band is performed by the software program (please contact Viro-Immun for further assistance). The result including the interpretation will be printed together with the protocol sheet.

**Note:**

Every [NTS] has a reaction control. If the assay handled right, you can see on top of the [NTS] the purple coloured reaction control (see control strip) .

The with of the bands as well as their intensity are variable and depend on the concentration of the detected Helicobacter pylori in the sample.

The band pattern in the IgG determination is generally more intense and diverse than IgA.

Helicobacter pylori reactive proteins

Table 1

Molekular weight in kDa	Antigen	Specificity
>125		not specific
p120	<b>Cag A</b> ( <i>Cytotoxin-associated-gene-A</i> )	<b>high specific</b> important pathogenicity factor for H. pylori, Ulcus marker
p117-91		not specific
p87	<b>VacA</b> ( <i>Vacuolating-Cytotoxin-A</i> )	<b>high specific</b> for H. pylori, frequent common with CagA, Ulcus marker
p84-68		not specific
p66	<b>UreB</b> ( <i>Big part of Urease</i> )	<b>specific for H. pylori</b> , marker for inflammation
p59	<b>HspB</b> ( <i>Heat shock protein B</i> )	<b>spezifisch</b> for H. pylori, cross reaction to other heat shock proteins known
p55	<b>FlaA</b> ( <i>Flagellin A-protein</i> )	<b>specific for H. pylori</b>
p50		<b>specific for H. pylori</b>
p37		<b>specific for H. pylori</b>
p33		<b>high specific</b> for H. pylori, functioning till now unknown
p29	<b>UreA</b> ( <i>small part of Urease</i> )	<b>specific for H. pylori</b> , marker for inflammation
p24		<b>high specific</b> marker for pathogenicity, functioning till now unknown

Helicobacter pylori IgG/IgA:

Both the intensity and position of bands should be taken into consideration when assessing WESTERN BLOT bands.

The **Helicobacter pylori** antigens **IgG** and **IgA** may be divided in WESTERN BLOT into three classes:

- A** **unspecific antigen bands**  
< 24kD / p68-p84kD / p91kD-p117kD / >125kD
- B** **specific antigen bands**  
UreA p29kD / p37kD / p50kD / FlaA p55kD/ HspB p59kD / UreB p66kD
- C** **highly specific antigen bands**  
**CagA** p120kD / **VacA** p87kD / p33kD / p24kD

**negative** - unspecific bands from A, weak antigen bands from B  
- 3 specific antigen bands from B  
- 1 highly specific band from C and 1 band from B

**questionable** - 4 specific antigen bands from B  
- 1 highly specific antigen band from class C and 2 specific bands from B  
- 2 clear highly specific bands from C

**positive** - 1 heavy and 1 further clearly highly specific band from C  
- at least 2 clear antigen bands from C and 1 specific from B  
- at least 1 clear band from C and at least 3 clear specific bands from B

unclear, weak bands should be not considered in evaluation.

**Diagnostic Relevance and interpretation of results**

- In the early phase of infection it could be happen that less antibodies will be detected, in this stage the test could be negative.
- If there is a clinical suspicion of Helicobacter infection and a negative or questionable test result, the tests should repeated after 2 weeks.

- A study shows that within the first weeks after treatment antibody titres had fallen by 20% – 30% irrespective of the success of bacterial eradication. In the bacteria-negative patients the decrease continued. 6 – 12 months after treatment the titre was 50% or less of pretreatment value. In patients who remained infected, the initial drops of antibody titres, if any, was followed by unchanged or slightly rising titres. IgA antibodies are valuable in patients who raise titres in this immunoglobulin class only.
- Seroconversion appears approximately 22 – 23 days after acute infection. For the interpretation of positive test results the clinical picture should always be taken into account. In case of a negative test result a florid H. pylori infection can almost be excluded. The majority of persons are infected lifelong but show no clinical symptoms.

IgG	IgA	Interpretation	Recommendation
-	-	no specific antibodies detectable,	For patients with increased risk of acute infection : monitoring examination.
+	+	acute or recent infection probable, In case of clinical symptoms treatment should be considered.	<b>Monitoring</b> of IgG, and IgA antibodies resp. use of other diagnostical parameters e.g. histological examination.
+- ++	-	<b>single class IgG response or endemic infection titre</b> In case of clinical symptoms treatment should be considered.	<b>Monitoring</b> of IgG, and IgA antibodies resp. use of other diagnostical parameters e.g. histological examination
-	+- ++	<b>single class IgA response possible</b> In case of clinical symptoms treatment should be considered.	<b>Monitoring</b> of IgG, and IgA antibodies resp. use of other diagnostical parameters e.g. histological examination

- Equivocal results** should be retested. Following the confirmation of the equivocal result the monitoring of the patient's antibodies is recommended in order to exclude unspecific reactions e.g. cross-reactivity, which may also cause equivocal results.

Note:

Reasons for **false negative results** are:

- low antibody response, 2. predominately development of IgA antibodies, 3. early stage of acute infection.

**False positive results** are based on cross reactions with e.g.: Campylobacter jejuni, E. coli, E. faecalis.

### Performance Characteristics

#### Specificity / Sensitivity

31 samples IgG and 37 samples IgA were tested parallel in Western Blot Anti-Helicobacter IgG and IgA and comparison methods ELISA. The sensitivity and specificity are based on the results found.

**Specificity:** IgG 100%      **Sensitivity:** IgG 90 %  
**Specificity:** IgA 100%      **Sensitivity:** IgA 94,1%

#### Cross-reactivity

- In the early infection phase there may be an insufficient number of antibodies and therefore a negative test result is possible.
- When there is a clinical suspicion of Helicobacter infection with a negative or questionable test result, a repeat test two weeks later is recommended.
- Antibiotic therapy may suppress the production of detectable antibodies in the early stages.
- The serological data should be taken into consideration with the clinical history for a full diagnosis. Where the results are questionable repeated testing is recommended.

#### Precision and reproducibility

Intra-assay reproducibility was determined by testing samples of different levels of antibody activities in one test run.

In determination of the samples there is no significant difference of detected antigen bands and their intensity.

Inter-assay reproducibility was determined by testing samples of different levels of antibody activities in different test runs.

In determination of the samples there is no significant difference of detected antigen bands and their intensity.

#### Internal quality control

The negative control controls should be carried out with every test run. Recommendation: For every new Lot, it is recommended to perform an internal quality control by testing one or more internal control samples which should be clearly specified.

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### Troubleshooting

ERROR	POSSIBLE CAUSES
no colourimetric reaction after addition substrate	no [CONJ]ALK pipetted, contamination of [CONJ]ALK (possibly with control sera during pipetting) may cause an inactivation.
generally too high reaction	incorrect [CONJ]ALK (i.e. not from original test kit), incubation time too long or incubation temperature too high, water quality for [WASHBUF] insufficient (low grade of deionization)
generally too low reaction	incorrect [CONJ]ALK (i.e. not from original test kit), incubation time too short, incubation temperature too low
false positive / negative samples	incorrect dilution of samples, lipaemic, haemolytic or icteric samples
unexplainable outliers	contamination of pipettes, tips or containers or with metals (iron, copper etc.),
high variation (from series to series)	incubation conditions not constant (time, temperature) high variation of room temperature, controls and samples are not carried out at same time (same intervals) check pipetting order, person related variation.
general source of error, that may leads to incorrect results	- surface damaging of [NTS] because of improper treatment - [NTS] dried out after washing (unreproducible results). - Reagents not pre-warmed to room temperature prior to use ([NTS])!. - Insufficient washing

The close observance of the test procedure will help to avoid these errors

### WESTERN BLOT SHORT ASSAY PROCEDURE

Dilute **patient sample 1:101** :  
10 µl patient sample + 1000µl [SPE]DIL.

Dilute **controls 1:10**:


900µl [SPE]DIL + 100 µl control

Dilute [CONJ]ALK[10x] 1:10:

f.e: 4,5ml [SPE]DIL + 0,5 ml [CONJ]ALK[10x]

Take required [NTS] holding with tweezers by the numbered end out of the tube.

All incubations / wash steps on the rocking platform

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