

## Aldosterone

Ref: AD-ALDO-(R1)

### Polyclonal rabbit Aldosterone-3CMO-BSA

Alphadia has others lots but cross-reactions are not yet established

#### Storage condition

- Buffer: phosphate 0.1M / pH 7.4
- Preservative: NaN3 1g/l
- Storage temperature: frozen (-20°C)

#### Cross reactivity were tested under next conditions

- Liquid phase: RIA
- Working dilution used: ND

#### Cross-reactivity

#### AD-ALDO-(R1)

3b, 5a Tetrahydroaldosterone	16.20%
3a, 5b Tetrahydroaldosterone	0.13%
Prednisolone	0.018%
Cortisol	<0.003%
11-Deoxycortisol	<0.003%
11-Deoxycorticosterone	<0.003%
17-OH Progesterone	<0.003%
Progesterone	<0.003%
Prednisone	<0.002%
Testosterone	<0.002%
Androstenedione	<0.002%