



**Progesterone Rabbit anti Human (AD-PROG3CMO)**

Lot: M82090376

Immunogen: 3CMO-BSA

Storage condition: -20°C

Cross reaction:

Liquid Phase: RIA

Working dilution: used 1/75.000

Sensitivity of the standard curve: 4 – 8 pg.

11 $\alpha$ -Hydroxyprogesterone	1%
17 $\alpha$ -Hydroxyprogesterone	0.08%
20 $\beta$ Hydroxyprogesterone	0.72%
20 $\alpha$ -Hydroxyprogesterone	0.06%
21 Hydroxyprogesterone	3.16%
17 $\beta$ oestradiol	less than 0.05%
Estradiol	less than 0.009%
Prednisolone	less than 0.05%
Pregnenolone	1.6%
Testosterone	0.08%
Dehydroepiandrosterone	0.08%
Cortisone	less than 0.05%
Dexamethasone	less than 0.05%
Cholesterol	less than 0.05%
Cortisol	less than 0.05%
11-Deoxycortisol	0.12%
21-Deoxycortisol	0.11%
11-Deoxycorticosterone	2.91%
Corticosterone	less than 0.05%
Androstenedione	0.17%

ALPHADIA S.A/N.V.

Siège d'exploitation : Rue du bosquet, 15A • BE-1435 Mont-Saint-Guibert • BELGIUM

Siège social : Rue René Descartes, 19 • BE-7000 Mons • BELGIUM

RPM Mons • TVA BE 0454.402.636

Tel : +32 (0) 10 68 56 10 • E-mail : [contact@alphadia.be](mailto:contact@alphadia.be) • [www.alphadia.be](http://www.alphadia.be)

DHEA  
DHEA-SO4

less than 0.05%  
less than 0.05%

ALPHADIA S.A/N.V.

Siège d'exploitation : Rue du bosquet, 15A • BE-1435 Mont-Saint-Guibert • BELGIUM

Siège social : Rue René Descartes, 19 • BE-7000 Mons • BELGIUM

RPM Mons • TVA BE 0454.402.636

Tel : +32 (0) 10 68 56 10 • E-mail : [contact@alphadia.be](mailto:contact@alphadia.be) • [www.alphadia.be](http://www.alphadia.be)